accomplished quickly and without delay in any step. Once the dried material has been reconstituted, the sample should be used immediately, as clotting will occur which could affect results.

Test procedure for HEMOCHRON F7CA510, FTK-ACT, APTT (A103), Citrate APTT (A104), PT (A201), or Citrate PT (202)
1. Peel back the flap of the crimp seal on the diluent vial ensuring that the level on the diluent vial matches the level on the Whole Blood CQ vial, but do not remove the entire crimp seal (use this technique when opening each vial).
2. Draw 20 cc of the crimp seal into the diluent vial using the syringe and needle. Transfer the diluent into the whole blood vial by direct puncture of the stopper.
3. Remove the syringe and needle from the whole blood vial. Do not discard syringe and needle.
4. Using the same syringe, withdraw 2.0 cc of the rehydrated whole blood control. Immediately transfer the 2.0 cc into the appropriate HEMOCHRON tube; the control material first enters the tube, depress the START button on the appropriate HEMOCHRON test well.
5. Immediately agitate the tube vigorously from end to end 10 times.
6. Insert the tube into the appropriate HEMOCHRON test well. Quickly rotate the tube clockwise until the detector green light is on. See the appropriate instrument operator manual for additional information.
8. Upon test completion, record results.

INTENDED USE
The HEMOCHRON® Whole Blood Quality Controls are intended to be used to perform quality control assays using the HEMOCHRON F7CA510, FTK-ACT, APTT, Citrate APTT, PT, Citrate PT, TT and HNTT test tubes.
For *in vitro* Diagnostic Use

SUMMARY AND EXPLANATION
Blood coagulation instruments and assays should be quality controlled prior to during routine use. A proper quality assurance program includes regular instrument maintenance, quality control assay measurements and complete performance records.
REAGENTS
The HEMOCHRON Whole Blood Quality Controls contains both Level 1 and Level 2 quality control materials. Assayed clotting times are provided with each lot of material.

Each Whole Blood Quality Control consists of:
• Whole blood controls, consisting of dried fixed bovine red blood cells, buffered sheep and horse plasma in individual test vials.
• Diluent preparation, consisting of distilled water, sodium chloride, Tween 20, calcium chloride, anticoagulant and preservatives in individual vials for reconstitution of the whole blood control.
NOTE: Diluents are not matched to the whole blood vials. DO NOT use diluents from one lot of control with whole blood vials from a different lot.
Level 1 and Level 2 quality control products are available for use with the following HEMOCHRON tests in the HEMOCHRON Model 401, 801 and Response instruments.
• Level 1 and Level 2 quality control, cat. # QCACT: for use with HEMOCHRON kaolin and Celite® ACT tubes.
• Level 1 and Level 2 quality control, cat. # QCAPT: for use with HEMOCHRON Citrate APTT tubes.
• Level 1 and Level 2 quality control, cat. # QCPT: for use with HEMOCHRON PT tubes.
• Level 1 and Level 2 quality control, cat. # QCCT: for use with HEMOCHRON TT and HNTT tubes.

CAUTION: Annual blood products are included in the Whole Blood Quality Control preparation. No material in any preparation is derived from human blood. All blood products, whether of human origin or animal origin, should be handled with care.

HANDLING AND STORAGE
When refrigerated (2°-8° C, 35-46°F), the vials are stable until the marked expiration date. The quality control product should never be exposed to temperatures in excess of 37° C (98.6°F), or below freezing. Reconstituted vials should be used immediately. The Whole Blood Quality Controls may also be stored at room temperature (15-30°C, 59-86°F) for up to 4 weeks but must never exceed the marked expiration date. If stored at room temperature, reading a 4 weeks on the enclosed label is necessary. The label should be affixed to the outside of the product box, covering the stamped expiration date.
PROCEDURE
Materials Provided
Whole Blood Quality Control kits for use with the HEMOCHRON F7CA510, FTK-ACT, APTT, Citrate APTT, PT, Citrate PT, TT and HNTT tubes contain the following:
• 10 vials of dried Level 1 whole blood control (2.0 cc).
• 10 vials of dried Level 2 whole blood control (2.0 cc).
• 10 vials of Level 2 diluent (4.0 cc).
• A quality control tracking chart.
Materials Required But Not Provided
• HEMOCHRON Model 401, 801 or Response instrument.
• HEMOCHRON Test tubes:
 F7CA510 (F7CA510, black top Celite activated ACT)
 FTK-ACT (kaolin activated ACT)
 APTT (A103 for fresh whole blood or A104 for citrated whole blood)
 PT (A201 for fresh whole blood or A202 for citrated whole blood)
 TT (A301)
• 3 cc plastic syringe with a 20 or 21 gauge needle for control vial rehydration.
Reagent Preparation
A. Quality Control Product
Remove the appropriate Whole Blood Quality Control vials from the refrigerator and allow them to come to room temperature (15-30°C, 59-86°F) prior to testing. This could require up to 60 minutes.
B. HEMOCHRON test tubes
Remove the appropriate test tubes from the refrigerator and allow them to come to room temperature prior to testing. This could require up to 60 minutes.

CAUTION: All used tubes and control vials should be considered as potentially infectious, handled with care, and disposed of in accordance with standard medical waste disposal policy.
NOTE: Reconstitution and mixing of the whole blood material should be performed immediately after reconstitution.

PERFORMANCE CHARACTERISTICS
Expected Values
The acceptable range for each Whole Blood Quality Control product is published on the enclosed label. These values have been established at the manufacturer's facility and are lot specific. The range represents the acceptable performance range for all HEMOCHRON tube instruments. The user should be able to obtain a value within the published range included for each product. A proper quality assurance program includes regular control measurements and maintenance of complete performance records.
ITC recommends that each institution establish its own expected range of response based on the mean ± 2 standard deviation of at least 20 repeated tests. The local mean values established should be compared to the manufacturer's acceptable performance range. Studies show that intra-laboratory variation in test results should produce a coefficient of variation of approximately 10% or less for coagulation control tests.
LIMITATIONS
Quality control material can significantly alter performance characteristics. Care must be taken to establish and adhere to the recommended test procedure found in this package insert and the associated test tube package insert. Variations in ambient temperature should be minimized for consistency of test performance. In cases where results fall outside of the assigned range, patient results should be considered suspect and performance of quality control tests should be repeated. Each institution should establish a tracking record to identify instruments, tubes, reagents or operators whose performance is outside of the institution's expected range.
In cases where results are outside of the expected range, the cause will likely fall into one of the following possible categories:
1. Operator technique
2. Quality Control Product/Test Tube
3. Instrumentation
4. Environmental factor

CAUTION: All used tubes and control vials should be considered as potentially infectious, handled with care, and disposed of in accordance with standard medical waste disposal policy.
NOTE: Reconstitution and mixing of the whole blood material should be performed immediately after reconstitution.

IF OUT-OF-RANGE RESULTS OR INSTRUMENT FAULT MESSAGES PERSIST, CONTACT ITC TECHNICAL SERVICES.
TROUBLESHOOTING
Please refer to the appropriate instrument Operator's Manual "Troubleshooting" section for instrument software fault messages, definitions and actions.

PROBLEM	CAUSE	ACTION
Control value is above or below published range	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.
	Time period between control material mixing addition to test tube is too long	Repeat test. After thorough mixing of the control material and diluent, IMMEDIATELY withdraw the rehydrated whole blood control and transfer to the appropriate HEMOCHRON tube.
Fault message was obtained (see appropriate instrument Operator's manual)	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

REFERENCES
1. Stalland, B and Westgard, J. Quality Control: Theory and Practice, "In Clinical Diagnosis and Management by Laboratory Methods", ed. Henry, J (Philadelphia WB Saunders Co., 1984) pp.74-93.
2. Colton, C. Statistics in Medicine (Boston: Little, Brown, and Co., 1974) pp 289-313.

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LITERATURANGABEN
1. Stalland, B and Westgard, J. Quality Control: Theory and Practice, "In Clinical Diagnosis and Management by Laboratory Methods", ed. Henry, J (Philadelphia WB Saunders Co., 1984) pp.74-93.
2. Colton, C. Statistics in Medicine (Boston: Little, Brown, and Co., 1974) pp 289-313.

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HEMOCHEON® Whole Blood Coagulation Systems

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