
# Hemochron®

## Whole Blood Coagulation System

### Protamine Dose Assay (PDAO)

## Package Insert

English

#### INTENDED USE

The Hemochron® PDAO is designed for in vitro diagnostic protamine dose determination. The PDAO is intended for use in the verification of heparin neutralization at the completion of cardiovascular surgery. The kaolin-activated PDAO assay may also be used for patients receiving aprotinin. The Hemochron PDAO should be used with established Hemochron procedural guidelines with Hemochron models 401, 801, 8000 and Response.

For In vitro Diagnostic Use / For Professional Use / Rx Only

#### SUMMARY AND EXPLANATION

The amount of protamine sulfate required to neutralize heparin after cardiovascular surgery is commonly based on a ratio of protamine to the amount of heparin infused during the procedure<sup>1</sup> or based on protamine titration.<sup>5,6</sup> Accurate determination of protamine dose is advantageous in view of the known complications associated with excess or insufficient protamine.<sup>6,8</sup>

Because of the pharmacokinetics of protamine, it is important to quantify and administer the correct protamine dosage required to neutralize the residual, circulating heparin.<sup>7</sup> Too low a dose may result in heparin-related bleeding and may require an additional protamine infusion to fully reverse the heparin effect.<sup>1,7</sup> Too high a dose of protamine, a more common occurrence, may result in bleeding due to the anticoagulant properties of protamine when not bound to the heparin molecule.<sup>1,7</sup> Infusion of excess protamine may also result in severe hypotension, anaphylactic shock, pulmonary vasoconstriction, and cardiac deceleration.<sup>9,10</sup>

The purpose of protamine titration is to determine the exact protamine dosage required to completely neutralize circulating heparin in any given patient. Protamine neutralization of heparin in vitro results in normalization of the clotting time in a manner which correlates with the in vivo effects of the drug.<sup>5,6</sup>

#### PROTAMINE DOSE ASSAY-O

The PDAO provides a convenient method for protamine titration using the principles of a in vitro dose-response heparin neutralization and verification of heparin neutralization.<sup>11</sup> Titration analysis performed using PDAO quantifies the milligrams of protamine per milliliter of patient blood necessary to neutralize the patient's residual heparin load and return the Activated Clotting Time (ACT) to the preheparinized baseline value.

The PDAO test is available with either Celite® (diatomaceous earth) or kaolin activators. The Celite-activated PDAO test is not intended for use with the protease inhibitor aprotinin (Trasylol®, Bayer Corporation), which may be administered to reduce post-operative bleeding, especially during cardiopulmonary bypass surgery, and can prolong the Celite-activated PDAO. Kaolin is unaffected by moderate doses of aprotinin. The kaolin-activated PDAO can be used as an alternative for protamine dose response evaluation in patients receiving aprotinin. The post-surgical patient should be monitored with the ACT and PDAO to determine the need to infuse additional protamine to compensate for latent heparin effect, which may possibly result in heparin rebound.<sup>10</sup>

The PDAO test is performed using Hemochron whole blood coagulation models 401, 801, 8000 or Response. The PDAO contains the reagents required to perform an ACT test plus an appropriate amount of protamine sulfate to identify the presence of residual heparin and determine an appropriate dose (if any) for neutralization.

In 1975, Ball et al.<sup>1</sup> published the dose-response curve method, based on the ACT test,<sup>10</sup> to individualize heparin and protamine dosages. In 1983, Dutton et al.<sup>2</sup> described a method of protamine titration using the Hemochron 400 system. In this method, varying amounts of liquid protamine are added to samples of the patient's heparinized blood, and the ACT is determined for each sample. ACT results plotted against in vitro protamine concentration constitute a dose-response heparin neutralization curve. The exact protamine concentration required to normalize the ACT to its level prior to heparinization may be calculated from such a curve.

The PDAO uses patient whole blood in conjunction with a status ACT to construct an individualized protamine titration curve. The protamine dose is calculated based on these clotting times and the patient's blood volume. This can be calculated automatically by the Hemochron Response (software version 2.0 or higher) or 8000, or manually using the worksheet provided.

#### REAGENTS

The Hemochron PDAO test tubes are intended for use with fresh whole blood. The test tubes contain a lyophilized preparation of protamine sulfate with added activator, stabilizers and buffers.

- PDAO contains protamine sulfate (0.01mg), diatomaceous earth, stabilizers and buffers.
- PDAOK contains protamine sulfate (0.01mg), kaolin, stabilizers and buffers.

#### PRODUCT USE WARNING

**NOTE: Observe universal precautions at all times.**

- The blood specimen should be transferred using an appropriate transfer needle to pierce the stopper.
- Always use a two-hand technique to transfer blood. One hand securely holds the tube while the second hand dispenses the blood specimen from the syringe.
- PDAO test tubes are made of glass. They can be broken or cracked if mishandled. Do not drop or toss tubes.
- The PDAO test tubes contain a particulate activator (Celite or kaolin) as well as a material of biological origin (protamine sulfate). Do not handle the material, aerosol or ingest.
- All used test tubes containing human derived blood should be discarded in approved biohazard containers.
- Do not remove stoppers to deliver a blood sample to the tube.

#### STORAGE & STABILITY

The Hemochron PDAO test tubes should be stored at controlled room temperature (15-30°C). The PDAO tubes are stable until the marked expiration date. Avoid prolonged exposure of the PDAO kits to temperatures exceeding 30°C.

#### SPECIMEN COLLECTION AND HANDLING

##### Materials Provided

- PDAO test tube (orange-top; Celite or kaolin-activated as appropriate).
- PDAO worksheet (for manual calculations, if required).

#### Materials Required

- Hemochron Response, 8000, 801 or two 401s.
- Non-siliconized, non-heparinized 5 mL syringes for blood collection.
- #20 gauge or larger straight needle for transfer of the blood specimen into the test tube.
- ACT test tube, HRF7CA510 (black top) or HRF7K-ACT (gold top), as appropriate.

**CAUTION: Every precaution should be taken to use proper technique with syringes to avoid a mLidental needlestick.**

For blood collection, adhere to the appropriate technique (A or B):

- In-dwelling venous or arterial blood line. (Do not obtain blood from a heparinized access line or indwelling heparin lock.)
    - Discontinue fluids drip, if required.
  - Using a two syringe technique, draw two 5 mL syringes and discard the first syringe.
- B. Extracorporeal blood line port.
- Flush the extracorporeal blood access line by withdrawing and discarding 5 mL of blood.
  - Draw 5 mL whole blood into a 5 mL syringe.

#### QUALITY CONTROL (QC)

Routine quality control testing and tracking should be a part of a comprehensive quality assurance program. Hemochron Whole Blood Coagulation System Quality Control products are available to make routine QC convenient and affordable.

#### Daily QC of the Instrument

Hemochron instruments should be quality controlled at two levels of control once every 8 hours of operation. To assist in accomplishing daily QC, Electronic System Verification Tubes are available and can provide multiple level (normal and abnormal) quality control checks on the instrument. A Temperature Verification Tube is also provided for use every 6 months to verify that the instrument is maintaining the proper temperature 37°C ± 1°C. Any errors found in the temperature controls system will also be displayed on the LCD screen by the instrument.

#### QC of Hemochron Test Tubes

Each box of Hemochron test tubes should be validated at least once prior to use. This can be accomplished by using the appropriate Hemochron Liquid Quality Control. Acceptable performance ranges for the test tubes are included with each Hemochron Quality Control Product Kit. After each individual box of test tubes has been verified, the "Performance Verified" table provided on the side panel of each test tube box should be completed. This box is now "IN CONTROL" and will not require further liquid quality control unless a shift in clinical results is experienced.

**NOTE: If multiple boxes are received within the same shipment, it is recommended to validate each box prior to use.**

#### TEST PROCEDURE

If using the Hemochron Response (software version 2.0 or higher) or 8000, the system software will quantify the protamine dose (if any) based on the PDAO test result and information provided by the user, along with the status ACT results. If using the Hemochron Response (software versions below 2.0), 401 or 801, use the PDAO calculation worksheet provided with the test kit to calculate the protamine dose.

If using the Celite-activated PDAO, an HRF7CA510 (black top) ACT should be used. If using the kaolin-activated PDAO, an HRF7K-ACT (gold top) should be used.

#### RECOMMENDED PROTOCOL

- Five minutes following initial protamine infusion, simultaneously perform a status ACT and a PDAO. From the collection syringe, dispense exactly 2.0 mL of blood into the PDAO test tube. At the same time, depress the **START** button of the appropriate Hemochron test well.
- Immediately agitate the test tube vigorously from end to end ten times.
- Insert the PDAO test tube into the appropriate test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- From the same collection syringe, dispense exactly 2.0 mL of blood into the ACT tube. At the same time, depress the **START** button of the appropriate Hemochron test well.
- Immediately agitate the test tube vigorously from end to end ten times.
- Insert the ACT test tube into the appropriate Hemochron test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- At the indicator tone, the test result.
- Determine the protamine infusion dose. See Dose Calculations below.

#### Dose Calculations

If using the Hemochron Response (software version 2.0 or higher), refer to the Hemochron Response RxDx Operator Manual for detailed instructions on the use of this automated calculation system. If using the Hemochron 8000, follow instructions in this instrument's operator manual. If using the Hemochron system with software versions below 2.0, Hemochron 401 or 801, follow instructions on the PDAO worksheet provided in this kit to calculate the protamine dose.

#### Verification of heparin neutralization

Residual heparin is reduced by an ACT result more than 10% higher than the concomitant PDAO result. A PDAO result which is higher than the corresponding ACT also indicates complete neutralization of circulating heparin.

#### LIMITATIONS OF THE PROCEDURE

The PDAO test is affected by poor technique in the test procedure. Proper specimen/reagent mixing is required for precise and accurate testing.

The accuracy of the protamine dose prediction is dependent upon the accuracy of the patient's calculated blood volume. In some instances, blood volume changes during surgery may be dramatic due to excessive intravenous-fluid infusion, urinary loss, blood loss, and blood component transfusion. As with any medication by infusion, the user must be aware of these alterations because of their potential effects on the efficacy of the protamine infusion.

The use of additional heparin or heparinized blood products (i.e., salvaged red blood cells from the extracorporeal system) following the performance of the PDAO will negate the test results and protamine prediction. In this circumstance, the PDAO assay should be repeated to reflect the additional protamine requirements.

Only kaolin-activated PDAO tubes should be used for patients receiving aprotinin. Studies have not been performed utilizing other antifibrinolytic agents.

The PDAO test is not intended for use when the status ACT exceeds 300 seconds.

6. Test results which do not agree with expected values should be verified and thereafter evaluated by alternative diagnostic methods.

#### PERFORMANCE CHARACTERISTICS

The PDAO was shown in laboratory and clinical studies to detect residual heparin above 0.25 units/mL and to quantify additional millinotes of protamine (as low as 25 mg) required to obtain complete heparin neutralization.<sup>12</sup>

#### PRECISION

Reproducibility was determined using Hemochron Quality Control PDAO.

PDAO	Level_1				Level_2			
	Day 1	Day 2	Day 3	Combined	Day 1	Day 2	Day 3	Combined
N	3	3	3	9	3	3	3	9
Mean	126	126	121	124	331	330	345	335
Std dev	3.5	3.1	3.2	3.6	0	4.0	13.2	10.0
CV%	2.7	2.4	2.6	2.9	0	1.2	3.8	3.0
PDAOK								
N	5	5	4	14	5	5	4	14
Mean	122	134	126	128	391	405	403	399
Std dev	9.7	3.9	10.6	9.4	13.3	13.9	15.7	14.6
CV%	7.9	2.9	8.4	7.4	3.4	3.4	3.9	3.7

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